

# Acetylcholine Assay Kit (Colorimetric/Fluorometric)

LS-K244-100 (100 Tests) • Store at -20°C



---

## Introduction

ACETYLCHOLINE is a neurotransmitter produced in acetylcholinergic neurons. It plays important roles in skeletal muscle movement, regulation of smooth and cardiac muscles, as well as in learning, memory and mood. This Method provides a simple, direct and high-throughput assay for measuring acetylcholine in biological samples. In this assay, acetylcholine is hydrolyzed by acetylcholinesterase to choline which is oxidized by choline oxidase to betaine and H<sub>2</sub>O<sub>2</sub>. The resulting H<sub>2</sub>O<sub>2</sub> reacts with a specific dye to form a pink colored product. The color intensity at 570nm or fluorescence intensity (530/585 nm) is directly proportional to the acetylcholine concentration in the sample.

## Key Features

- Use 20 µL samples. Linear detection range: colorimetric assay 10 to 200 µM, fluorometric assay 0.4 to 10 µM acetylcholine.

## Applications

- Assays: acetylcholine in biological samples such as serum, plasma, urine, saliva, milk, tissue, and cell culture.
- Drug Discovery/Pharmacology: effects of drugs on acetylcholine metabolism.

## Components

Component	K244-100
	100 Tests
Assay Buffer	10 mL
Enzyme Mix	Dried
ACHE Enzyme	120 µL
Dye Reagent	120 µL
Standard (2 mM acetylcholine)	400 µL

## Materials Not Supplied

Pipetting devices, centrifuge tubes, clear flat-bottom uncoated 96-well plates, optical density plate reader; black flat-bottom uncoated 96-well plates, fluorescence plate reader.

## Storage

The kit is shipped on ice. Store all components at -20°C. Shelf life of six months after receipt.

---

**FOR RESEARCH USE ONLY! Not for use in humans.**

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121  
[www.LSBio.com](http://www.LSBio.com) • (206) 464-1554 • [TechnicalSupport@LSBio.com](mailto:TechnicalSupport@LSBio.com)

# Acetylcholine Assay Kit (Colorimetric/Fluorometric)

LS-K244-100 (100 Tests) • Store at -20°C



## Assay Procedure

### Colorimetric Procedure

#### Sample treatment

Liquid samples such as serum and plasma can be assayed directly. Tissue and cell lysates can be prepared by homogenization in cold 1 x PBS and centrifugation (5 min at 14,000 rpm). Use clear supernatants for assay. Milk samples should be cleared by mixing 600  $\mu$ L milk with 100  $\mu$ L 6 N HCl. Centrifuge 5 min at 14,000 rpm. Transfer 300  $\mu$ L supernatant into a clean tube and neutralize with 50  $\mu$ L 6 N NaOH. The neutralized supernatant is ready for assay (dilution factor  $n = 1.36$ ).

Note: (1). SH-containing reagents (e.g.  $\beta$ -mercaptoethanol, dithiothreitol,  $> 5 \mu$ M) are known to interfere in this assay and should be avoided in sample preparation. (2). This assay is based on an enzymecatalyzed kinetic reaction. Addition of Working Reagent should be quick and mixing should be brief but thorough.

#### Procedure

1. Equilibrate all components to room temperature. Briefly centrifuge the tubes before opening. Keep thawed tubes on ice during assay. Reconstitute Enzyme Mix with 120  $\mu$ L Assay Buffer. Reconstituted Enzyme Mix is stable for 1 month when stored at -20°C. Note: a yellow precipitate may form after thawing reconstituted Enzyme Mix. If a precipitate forms, pellet it by centrifuging for 2 min at 14000 rpm and use the clear supernatant.
2. Standards: mix 24  $\mu$ L 2 mM Standard with 216  $\mu$ L dH<sub>2</sub>O (final 200  $\mu$ M). Dilute standard in dH<sub>2</sub>O as follows.

No	200 $\mu$ M STD + H <sub>2</sub> O	Vol ( $\mu$ L)	Acetylcholine ( $\mu$ M)
1	100 $\mu$ L + 0 $\mu$ L	100	200
2	60 $\mu$ L + 40 $\mu$ L	100	120
3	30 $\mu$ L + 70 $\mu$ L	100	60
4	0 $\mu$ L + 100 $\mu$ L	100	0

Transfer 20  $\mu$ L diluted standards into separate wells of a clear flat-bottom 96-well plate.

Samples: transfer 20  $\mu$ L of each sample into separate wells of the plate. Note: if a sample is known to contain choline, prepare an extra sample blank well with 20  $\mu$ L of the sample.

3. Color reaction. Prepare enough Working Reagent by mixing, for each well, 85  $\mu$ L Assay Buffer, 1  $\mu$ L ACHE Enzyme, 1  $\mu$ L Enzyme Mix and 1  $\mu$ L Dye Reagent. Add 80  $\mu$ L Working Reagent to each well. Note: for samples that contain choline, prepare a blank control reagent with no ACHE Enzyme (i.e., 85  $\mu$ L Assay Buffer, 1  $\mu$ L Enzyme Mix and 1  $\mu$ L Dye Reagent). Add 80  $\mu$ L of the control Reagent to each Sample Blank well.

Immediately tap plate to mix. Incubate 30 min at room temperature.

4. Read optical density at 570 nm (550-585 nm).

### Fluorometric Procedure

The fluorometric assay procedure is similar to the colorimetric procedure except that (1) 0, 3, 6 and 10  $\mu$ M acetylcholine standards and (2) a black 96-well plate are used. Read fluorescence intensity at  $\lambda_{ex} = 530$  nm and  $\lambda_{em} = 585$  nm.

Note: if the calculated acetylcholine concentration of a sample is higher than 200  $\mu$ M in the Colorimetric Assay or 10  $\mu$ M in the Fluorometric Assay, dilute sample in water and repeat the assay. Multiply result by the dilution factor  $n$ .

**FOR RESEARCH USE ONLY! Not for use in humans.**

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121

[www.LSBio.com](http://www.LSBio.com) • (206) 464-1554 • [TechnicalSupport@LSBio.com](mailto:TechnicalSupport@LSBio.com)

# Acetylcholine Assay Kit (Colorimetric/Fluorometric)

LS-K244-100 (100 Tests) • Store at -20°C

## Calculations

Subtract blank value (#4) from the standard values and plot the  $\Delta OD$  or  $\Delta F$  against standard concentrations. Determine the slope and calculate the acetylcholine concentration of Sample,


$$[\text{Acetylcholine}] = \frac{R_{\text{SAMPLE}} - R_{\text{BLANK}}}{\text{Slope } (\mu\text{M}^{-1})} \times n \quad (\mu\text{M})$$

$R_{\text{SAMPLE}}$  and  $R_{\text{BLANK}}$  are optical density or fluorescence intensity readings of the Sample and H<sub>2</sub>O Blank (or Sample Blank if sample contains choline), respectively.  $n$  is the sample dilution factor.

Conversions: 1 mM acetylcholine equals 14.6 mg/dL, 0.015% or 146 ppm.

## Sample Data

Acetylcholine Standard Curves



Version: V.08.09.2018

**FOR RESEARCH USE ONLY! Not for use in humans.**

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121

[www.LSBio.com](http://www.LSBio.com) • (206) 464-1554 • [TechnicalSupport@LSBio.com](mailto:TechnicalSupport@LSBio.com)